Fig. 1

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Electron Beam Diffraction Tube

The electron beam diffraction tube is used to make visible the interferences of electron beams on polycrystalline graphite and thus to demonstrate the wave properties of electrons.

Examples of Experiments:

- Diffraction at a polycrystalline lattice (Debye-Scherrer diffraction) *)
- Determination of wavelength
- Verification of de Broglie's equation $\lambda = \frac{h}{m \cdot v}$
- Structure investigations (interplanar spacings of graphite)

1 Safety Notes

- Protect the high-vacuum tube against fall, shock and similar. Risk of implosion!
- Check emission current I_{EM} and graphite foil during the experiment. If I_{EM} is higher than $200 \mu A$ or the foil glows up, immediately interrupt the connection to the anode voltage source (see Section 3).
- Cover plug pin ④ by safe-to-touch plug adapter (from 501 641 and following numbers).

2 Description, Standard Equipment, Technical Data

- ① Electrode system with indirectly heated cathode
- ② Carrier net, covered with graphite powder (polycrystalline diffraction grating)
- ③ Fluorescent screen
- ④ 4-mm plug pin, connected with the anode
- ⑤ Pair of 4-mm sockets, connected with heater filament
- ⑥ 2-mm socket connected with the cathode via internal resistor
- ⑦ Adapter 2-mm plug/4-mm socket for socket ⑥

Not shown:

Magnet for readjusting the direction of the electron beam should it fall on a spot on the carrier net ② not fully covered with graphite due to the manufacturing process or burn out.

*) Recommended for demonstrating Debye-Scherrer diffraction using visible light: Optical analogue to Debye-Scherrer interference (555 19).

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Heizspannung U_F : max. 6 V~
 Heizstrom I_F : max. 0,3 A
 Anodenspannung U_A : max. 5 kV~
 Emissionsstrom I_{EM} : max. 200 μ A

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
$$d_1 = 2.13 \cdot 10^{-10} \text{ m}$$

$$d_2 = 1.23 \cdot 10^{-10} \text{ m}$$

Distance between graphite grating and fluorescent screen: approx. 13.5 cm

Diameter of glass bulb: approx. 13 cm

3 Experiment Description

Debye-Scherrer diffraction of electron beams (determination of wavelength; verification of λ_{ex} by means of de Broglie's equation; determination of h; investigation of structure)

From the electrons emitted by a hot cathode a small beam is singled out through a pin diaphragm. After passing through a focusing electron-optical system, it is incident as sharply limited monochromatic beam on a polycrystalline graphite foil the atoms of which are arranged similar to a space lattice and act as a diffraction grating. The diffraction pattern as it appears on a fluorescent screen shows two concentric rings centered around the undiffracted electron beam. The diameter of the concentric rings changes with the wavelength, determined by the acceleration voltage, whereby each of the two rings corresponds to a Bragg reflection at the atoms forming a plane of the graphite.

For interference reflection of electron beams at lattice atoms, Bragg's equation

$$2 d \sin \vartheta = n \lambda \quad (1)$$

applies as for other waves too (e. g. X-rays).

(d = interplanar spacing of space lattice
 ϑ = reflection angle
 λ = wavelength of electron beams)

Hence, interference reflection occurs when the path difference $2 d \sin \vartheta$ equals an integral multiple of the wavelength.

If now electron beams strike the graphite foil, there will always be some angles with the incident radiation that fit Bragg's equation (1) due to the arbitrary arrangement of the microcrystals.

All reflections from such crystallites lie on cone-shaped shells with a common axis. The fluorescent screen of the tube represents a section through these cones so that the diffracted electron beams appear as concentric circles round the non-diffracted electron beam.


From the schematic representation of a reflected beam (Fig. 2) it follows that

$$\tan 2 \vartheta = r/l$$

From this it follows for small angles

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When substituting (2) into (1) it follows for the wavelength λ (with $n = 1$)

$$\lambda = \frac{r \cdot d}{l} \quad (3)$$

Supposing the wavelength λ , which can be determined from the acceleration voltage (see part 3 of the evaluation), is known, it is then possible to determine the interplanar spacings of graphite from the radii of the two diffraction rings:

$$d = \frac{\lambda \cdot l}{r} \quad (4)$$

Apparatus:	Cat. No.
1 Electron beam diffraction tube	555 17
1 Stand for electron tubes	555 05
1 High-voltage power supply unit, 10 kV	522 37
1 Filament voltage source, 6 V a.c., with current limiting to 0.3 A, e. g. Transformer 6 V/12 V, 30 W	562 73
Rheostat, 11 Ω	537 29
Measuring instruments for 300 μ A d.c., 0.3 A, e. g. Moving-coil instrument D	531 781
Demonstration multimeter	531 91
Double-pole cut-out switch	504 50
Plug adapter	from 501 641
Connecting leads	from 501 20 and following numbers

Experimental set-up (see Fig. 3):

Arrange the power supply units so that their transformers do not exert any magnetic influence on the electron beams; keep them at a distance.

It is advisable to carry out the experiment in a darkened room.

Carrying out the experiment:

Switch on filament voltage and limit the heater current I_F to 0.3 A.

Gradually increase the anode voltage from 0 to 5 kV max.

Important: Check emission current I_{EM} (max. 200 μ A) and graphite foil *continuously*.

If I_{EM} is higher than 200 μ A or the foil glows up, *immediately* interrupt the connection to the voltage source via the double-pole cut-out switch, reset U_A to 0 V and discharge power supply unit by shorting. Position the adjusting magnet on the tube neck (e. g. fix it with adhesive tape) so that when applying the anode voltage again (increase gradually) the electron beams fall onto another spot on the foil.

If the diffraction rings are not satisfactory, change the direction of the electron beams by means of the adjusting magnet so that they are incident on a spot of the graphite grating which is free of defects and two rings (1st order diffraction on grate planes with $d_1 = 2.13 \cdot 10^{-10}$ m and $d_2 = 1.23 \cdot 10^{-10}$ m) are clearly visible.

Determine diameters $2r_1$ and $2r_2$ of the diffraction rings at $U_A = 5$ kV (sliding caliper).

Determine by estimation the distance l (see Fig. 2) from the glass-bulb diameter of 13 cm.

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Evaluation:

1. Determination of wavelength

From the measured values for r_1, r_2 and l and the interplanar spacings of graphite

$$d_1 = 2,13 \cdot 10^{-10} \text{ m}$$
$$d_2 = 1,23 \cdot 10^{-10} \text{ m}$$

the resulting wavelength for an acceleration voltage $U_A = 5 \text{ kV}$ is

$$\lambda_1 = 0,176 \cdot 10^{-10} \text{ m and } \lambda_2 = 0,177 \cdot 10^{-10} \text{ m}$$

2. Comparison of λ_{ex} with the theoretical value resulting from the Broglie's equation

To verify the empirical value of λ , de Broglie's equation is applied:

$$\lambda = \frac{h}{m \cdot v}$$

When substituting the velocity of electrons v resulting from the equation of energy

$$\frac{1}{2} m v^2 = e \cdot U$$

for the electrons accelerated by a voltage U , then it follows:

$$\lambda = \frac{h}{\sqrt{2 e m U}} \quad (5)$$

$$= \frac{6,6256 \cdot 10^{-34} \text{ Js}}{\sqrt{2 \cdot 1,6021 \cdot 10^{-19} \text{ C} \cdot 9,1091 \cdot 10^{-31} \text{ kg} \cdot 5 \cdot 10^3 \text{ V}}}$$
$$= 0,173 \cdot 10^{-10} \text{ m}$$

The values of $0,176 \cdot 10^{-10}$ and $0,177 \cdot 10^{-10}$ as determined from the diffraction pattern agree quite well with the theoretical values.

3. Determination of h from the de Broglie equation

From (5) it follows that

$$h = \lambda_{\text{ex}} \sqrt{2 e m U} \text{ with } \lambda_{\text{ex}} = 0,1765 \cdot 10^{-10} \text{ m}$$

$$= 0,1765 \cdot 10^{-10} \text{ m} \cdot \sqrt{2 \cdot 1,6021 \cdot 10^{-19} \text{ C} \cdot 5 \cdot 10^3 \text{ V} \cdot 9,1091 \cdot 10^{-31} \text{ kg}} = 6,743 \cdot 10^{-34} \text{ Js}$$

(The comma has to be read as decimal point)

4. Determination of the interplanar spacings of graphite

From the measured values for r_1, r_2 and l as well as the wavelength λ which is calculated after eq. (4) from the acceleration voltage U_A

$$\lambda = \frac{h}{2 e m U_A} \quad \text{with } h = 6,6256 \cdot 10^{-34} \text{ Js}$$
$$e = 1,6021 \cdot 10^{-19} \text{ C}$$
$$m = 9,1091 \cdot 10^{-31} \text{ kg}$$
$$= 0,173 \cdot 10^{-10} \text{ m}$$

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